

# A high-throughput gut-immune co-culture model for identifying anti-inflammatory therapies

s. Peddibhotla, L. Boone, E. Taylor, B. McQueen, N. Murr, M. Bunger, E. Boazak. Altis Biosystems, Inc., Durham, NC


## Purpose and Objectives

Current treatments for inflammatory bowel disease (IBD) are limited in their long-term effectiveness, as many patients become non-responsive to anti-inflammatory drugs after several months of treatment, underscoring the urgent need for new therapeutic options. Animal models are poor tools for screening anti-inflammatory leads due to their low throughput, interspecies variability in immune response, and differences in gut physiology.

### Objectives



- Develop a human-relevant and high-throughput indirect gut-immune co-culture system that accurately mimics the human IBD microenvironment
- Demonstrate co-culture model responsiveness to endotoxin and prevention of damage with anti-inflammatory treatments

## Methods



- RepliGut® Planar – Transverse Colon was established on Transwell® membrane insets in 96-well format with standard methods
- THP-1 monocytes were directed to the macrophage lineage using PMA on bottom surface of standard 96-well plate.
- The two culture systems were brought together in the presence or absence of macrophage activation stimuli (IFN-γ and LPS) to create comparable “healthy” and “inflamed” intestine environments.


## Experimental Overview




- (A) Schematic description of gut-immune co-culture and treatment and assay timelines.
- (B) RepliGut® mean TEER values remained steady at 48h of co-culture with increasing THP-1 cell counts.

## Establishing Inflamed Conditions


### (A) CellTiter-Glo®



### (C) Cytokine Analysis



### (B) TEER




n=3 technical replicates. \*\*\*\*p<0.0001, Two-way ANOVA with Tukey's multiple comparison test


- RepliGut®-TC + THP1 co-cultures maintained high TEER and viability similar to RepliGut®-TC monoculture in the absence of IFN-γ and LPS.
- Pretreatment of THP-1 and co-treatment of RepliGut®-TC + THP-1 with IFN-γ and LPS caused a dose dependent decrease in (A) Cell Viability and (B) TEER after 48h of co-culture
- (C) Dose dependent increases in released cytokines with increasing pre-treatment of cultures by IFN-γ

## Preventing Inflammation

### (A) TEER



### (B) Cell Titer Glo®



n=3 technical replicates. \*\*\*\*p<0.0001, Two-way ANOVA with Tukey's multiple comparison test.

- TEER (A) and CellTiter-Glo (B) assays on inflamed tissues treated with Adalimumab (Ada, TNF-α inhibitor) and Tofacitinib (Tofa, Jak-Stat inhibitor)
- Ada and Tofa treatments fully prevented TEER reduction and cell viability loss

## Conclusions

- Human intestine-immune co-culture system established by co-culturing RepliGut® TC with THP-1-derived macrophages.
- Activation of macrophages to inflamed state using IFN-γ and LPS damaged RepliGut®-TC monolayer in the co-culture system and caused increased inflammatory cytokine release.
- Induction of inflammation in the RepliGut®-TC + THP-1 co-culture provides a tunable environment for comparing treatment efficacies in healthy and inflamed models.

## References

Pike, Colleen M., et al. "Characterization and Optimization of Variability in a Human Colonic Epithelium Culture Model". ALTEx - Alternatives to Animal Experimentation, vol. 41, no. 3, July 2024, pp. 425-38, doi:10.14573/altex.2309221.

Kämpfer AAM, Urban P, Gloria S, Kanase N, Stone V, Kinsner-Ovakainen A. Development of an in vitro co-culture model to mimic the human intestine in healthy and diseased state. Toxicol In Vitro. 2017 Dec;45(Pt 1):31-43. doi: 10.1016/j.tiv.2017.08.011. Epub 2017 Aug 12. PMID: 28807632; PMCID: PMC5744654.